FOOD SYSTEM

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FOOD SYSTEM

ENGLISH

**System for detection and presumptive identification
of pathogenic microorganisms from foodstuffs validated to ISO 16140 standard.**

DESCRIPTION

FOOD-SYSTEM is a 24-well system containing desiccated biochemical substrates and culture media for detection and presumptive identification of microorganisms from meat, milk and cheese, fish and other food products.

The system provides detection and presumptive identification of: *Salmonella* spp., *Proteus* spp. / *Providencia* spp., *Pseudomonas* spp., *Staphylococcus aureus*, *E. coli*, *Bacillus cereus*, *Listeria* spp., yeasts and moulds, and particularly it is validated to ISO 16140 standard for the detection of *Salmonella* spp. and *Listeria* spp.

The system is inoculated with a suspension of food sample and incubated at 36 ± 1 °C for 18-24 hours.

The tests for detection and presumptive identification of the microorganisms present in the sample are interpreted by assessing the colour change in the various wells.

KIT CONTENTS

The kit contains:

| | |
|---|---------------------------|
| 20 FOOD SYSTEM (40 tests) | 1 <i>TEST RESULT FORM</i> |
| 40 vials of <i>Physiological solution</i> (4.5 mL/vial) | 20 Sealing films |
| 1 Instruction sheet | |

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

| | |
|--|---|
| VASELIN OIL (ref. 80278) | H ₂ O ₂ Reagent (ref. 80057) |
| OXIDASE TEST STICK (ref. 88029) | BUFFERED PEPTONE WATER (ref. 20099, 412090, 611014) |
| KOVAC'S – REAGENT (ref. 80270) | RINGER'S <i>Solution</i> (ref. 81059) |
| Sundry materials used in microbiology laboratories; cover slides and slides; microscope. | |

CONFIGURATION

The configuration of the system is shown in Table 1. Each system allows two samples to be examined or one sample in double test.

Table 1:

| | |
|-------------------------|--|
| Well | PRESUMPTIVE IDENTIFICATION of <i>Salmonella</i> spp. |
| 1-LDC | Decarboxylation of lysine |
| 2-H₂S | Production of hydrogen sulphide |
| Well | PRESUMPTIVE IDENTIFICATION of <i>Proteus</i> spp. |
| 3-UR | Hydrolysis of urea |
| 4-PRO | <i>Proteus</i> spp. / <i>Providencia</i> spp. |
| Well | PRESUMPTIVE IDENTIFICATION of <i>Pseudomonas</i> spp. |
| 5-PSE | <i>Pseudomonas</i> spp. |
| Well | PRESUMPTIVE IDENTIFICATION of <i>Staphylococcus aureus</i> |
| 6-STA | <i>Staphylococcus aureus</i> |
| Well | PRESUMPTIVE IDENTIFICATION of <i>Escherichia coli</i> and enteropathogenic <i>E.coli</i> O157 |
| 7-ESC | <i>Escherichia coli</i> |
| 8-IND | Indole test for confirmation of <i>Escherichia coli</i> |
| Well | PRESUMPTIVE IDENTIFICATION of <i>Bacillus cereus</i> |
| 9-BCE | <i>Bacillus cereus</i> |
| Well | PRESUMPTIVE IDENTIFICATION of <i>Listeria</i> spp. |
| 10-LIS | <i>Listeria</i> spp. |
| 11-CAT | Catalase test for confirmation of <i>Listeria</i> spp. |
| Well | PRESUMPTIVE IDENTIFICATION of Yeasts and moulds |
| 12-Y/M | Yeasts and moulds |

PRINCIPLE OF THE METHOD

FOOD SYSTEM allows the detection and presumptive identification of the following microorganisms: *Salmonella* spp., *Proteus* spp. / *Providencia* spp., *Pseudomonas* spp., *Staphylococcus aureus*, *E. coli*, *E. coli* 0157, *Bacillus cereus*, *Listeria* spp., yeasts and moulds.

- The presence of ***Salmonella* spp.** is shown by a yellow to red-orange colour change of well **1-LDC** and by the yellow to black color change of well **2-H₂S**.

The confirmation of *Salmonella* spp. is performed by taking a drop of the culture broth from the well **2-H₂S** and plating it onto a selective agar medium (i.e. XLD agar, ref. 10056).

Proceed to identification of bacterial colonies according to ISO 6579:2004 ².

The yellow to black color change of well **2-H₂S** and no color change in the well **1-LDC** can indicate the presence of ***Citrobacter* spp.**

- The presence of ***Proteus* spp. / *Providencia* spp.** is shown by the yellow to red-fuchsia color change of the well **3-UR**, and by the yellow to brown-black color change of the well **4-PRO**.

- The presence of ***Pseudomonas* spp.** is shown by the yellow to turbid green colour change of the well **5-PSE**.

The confirmation of *Pseudomonas* spp. is performed by the oxidase test directly in the well (OXIDASE TEST STICK ref. 88029).

- The presence of ***Staphylococcus aureus*** is shown by the appearance of a black ring on the bottom of the well **6-STA**.

The confirmation of *Staphylococcus aureus* is performed by taking a drop of the culture broth from the well **6-STA** by means of a sterile loop, and plating it onto a selective agar medium (i.e. Baird Parker Agar + RPF, ref. 10521).

Proceed to identification of bacterial colonies according to ISO 6888-2:1999 ³.

- The presence of ***E. coli*** and ***E. coli* 0157** is shown by the red to blue color change of the well **7-ESC** and by the formation of a pink-red coloring following the addition of Kovac's reagent into the well **8-IND**.

- The presence of ***Bacillus cereus*** is shown by the yellow to turbid green color change of the well **9-BCE**.

The confirmation of *Bacillus cereus* is performed by taking a drop of the culture broth from the well **9-BCE** by means of a sterile loop, and plating it onto a selective agar medium (i.e. Bacillus Cereus Agar-MYP, ref. 10027).

Proceed to identification of bacterial colonies according to ISO 7932:2004 ⁴.

- The presence of ***Listeria* spp.** is shown by the yellow to black color change of the well **10-LIS** and by the development of bubbles following the addition of Reagent H₂O₂ into the well **11-CAT**.

The confirmation of *Listeria* spp. is performed by taking a drop of the culture broth from the well **10-LIS** by means of a sterile loop, and plating it onto a selective agar medium (i.e. O.A. Listeria Agar, ref. 10620).

Proceed to identification of bacterial colonies according to ISO 11290:2002 ⁵.

- The presence of **yeasts and moulds** is shown by the green to yellow color change of the well **12-Y/M** and by watching for chlamydospores and hyphae at the microscope.

COMPOSITION

Table 2:

| Well | Contents |
|--------------------|--|
| 1-LDC | Culture medium with substrate for the decarboxylation of lysine |
| 2-H ₂ S | Culture medium with substrate for the production of hydrogen sulphide |
| 3-UR | Culture medium with substrate for hydrolysis of urea |
| 4-PRO | Culture medium with substrate for the isolation of <i>Proteus/Providencia</i> spp. |
| 5-PSE | Culture medium with substrate for the isolation of <i>Pseudomonas</i> spp. |
| 6-STA | Culture medium with substrate for the isolation of <i>Staphylococcus aureus</i> |
| 7-ESC | Culture medium with substrate for the isolation of <i>Escherichia coli</i> |
| 8-IND | Culture medium to show the production of bubbles |
| 9-BCE | Culture medium with substrate for the isolation of <i>Bacillus cereus</i> |
| 10-LIS | Culture medium with substrate for the isolation of <i>Listeria</i> spp. |
| 11-CAT | Culture medium to show the catalase reaction |
| 12-Y/M | Culture medium with substrate for the isolation of yeasts and moulds |

Physiological solution (g/l): Sodium chloride 9 g, Distilled water 1000.0 mL, pH 7.0 ± 0.2

COLLECTION AND STORAGE OF THE SAMPLES

The samples to be subjected to the test must be collected in accordance with the standard methods for samples to be subjected to microbiological examination and must be sent to the laboratory for seeding into the **FOOD SYSTEM** immediately after being taken.

TEST PROCEDURE**a) PREPARATION OF THE SAMPLE**

1. Homogenize a quantity of foodstuff (1, 10 or 25 g) in a proper volume (9, 90 or 225 mL) of BUFFERED PEPTONE WATER or RINGER'S Solution in accordance with the standard procedures for the food sample in question. The final dilution of the sample should be 1:10. Store the homogenized sample at 2-4 °C until test ending.
2. Take 10 mL of homogenized sample (dil. 1:10) and transfer into a suitable tube.
3. Cover the tube and incubate at 36 ± 1 °C for 12-18 hours.

b) INOCULATION OF THE SYSTEM

1. Take a system from its wrapper and allow it to come to room temperature.
2. Write down the sample data and the starting date of the examination.
3. Transfer 0.5 mL, from the pre-incubated tube containing sample suspension, into a vial of *Physiological Solution* contained in the kit.
4. Distribute 0.2 mL (4-5 drops) of suspension in each well of the system.
5. Cover the first 3 wells **1-LDC**, **2-H₂S** and **3-UR** with 2 drops each of Vaseline oil.
6. Cover the system with the lid provided and incubate at 36 ± 1 °C for 18-24 hours.

The system is configured for performance of two tests. In case of performance of one test, cover the unused half system with the sealing film provided in the kit.

INTERPRETATION OF THE RESULTS

At the end of the incubation:

- Add 2 drops of Kovac's Reagent to the well **8-IND** and watch for the appearance of a pink-red ring.
- Add 2 drops of Reagent H₂O₂ to the well **11-CAT** and watch for the formation of bubbles.
- Observe the color changes in of the wells and interpret the results according to the Table 3.
- Note the results on the *TEST RESULT FORM* (copy as many forms as necessary).
-

Table 3:

| Well | Presumptive Identification of <i>Salmonella</i> spp. | Well color | |
|-------------------------|--|---------------------------------------|--------------------------------------|
| | | Positive reaction | Negative reaction |
| 1-LDC | Decarboxylation of lysine | red-orange | yellow |
| 2-H₂S | Production of hydrogen sulphide | black | yellow |
| Well | Presumptive Identification of <i>Proteus</i> spp. / <i>Providencia</i> spp. | Well color | |
| | | Positive reaction | Negative reaction |
| 3-UR | Hydrolysis of urea | red-fuchsia | yellow |
| 4-PRO | <i>Proteus/ Providencia</i> spp. | brown-black | yellow |
| Well | Presumptive Identification of <i>Pseudomonas</i> spp. | Well color | |
| | | Positive reaction | Negative reaction |
| 5-PSE | <i>Pseudomonas</i> spp. | turbid green | yellow |
| Well | Presumptive Identification of <i>Staphylococcus aureus</i> | Well color | |
| | | Positive reaction | Negative reaction |
| 6-STA | <i>Staphylococcus aureus</i> | black | yellow |
| Well | Presumptive Identification of <i>Escherichia coli</i> , <i>Escherichia coli</i> 0157 | Well color | |
| | | Positive reaction | Negative reaction |
| 7-ESC | <i>Escherichia coli</i> , <i>Escherichia coli</i> 0157 | blue | grey-red |
| 8-IND | Indole test | pink-red | yellow |
| Well | Presumptive Identification of <i>Bacillus cereus</i> | Well color | |
| | | Positive reaction | Negative reaction |
| 9-BCE | <i>Bacillus cereus</i> | turbid green | yellow |
| Well | Presumptive Identification of <i>Listeria</i> spp. | Well color | |
| | | Positive reaction | Negative reaction |
| 10-LIS | <i>Listeria</i> spp. | black | yellow |
| 11-CAT | Catalase test | presence of bubbles | absence of bubbles |
| Well | Presumptive Identification of Yeasts and Moulds | Well color | |
| | | Positive reaction | Negative reaction |
| 12-Y/M | Yeasts and Moulds | yellow | verde |
| | Microscopic observation | presence of chlamydospores and hyphae | absence of chlamydospores and hyphae |

QUALITY CONTROL

Each batch of **FOOD SYSTEM** is subjected to quality control using the following reference strains of microorganisms:

| | | | | | |
|-------------------------------|------------|-------------------------------|------------|---------------------------------|------------|
| <i>Salmonella typhimurium</i> | ATCC 14028 | <i>Pseudomonas aeruginosa</i> | ATCC 27853 | <i>Bacillus cereus</i> | ATCC 11778 |
| <i>Citrobacter freundii</i> | ATCC 8090 | <i>Staphylococcus aureus</i> | ATCC 25923 | <i>Listeria monocytogenes</i> | ATCC 35152 |
| <i>Proteus mirabilis</i> | ATCC 25933 | <i>Escherichia coli</i> 0157 | ATCC 35150 | <i>Saccharomyces cerevisiae</i> | ATCC 9763 |

PERFORMANCE

The validation of **FOOD SYSTEM** has been represented by a comparative study with the procedure indicated in the ISO 11290 for the isolation of *Listeria monocytogenes*, and with the procedure indicated in the ISO 6579 for the isolation of *Salmonella* spp. Overall 712 samples of meat and dairy products have been tested in duplicate with **FOOD SYSTEM** and reference method. The obtained results are indicated in the following table:

| | <i>Listeria</i> spp. detection | | <i>Salmonella</i> spp. detection | |
|----------------------|--------------------------------|----------------|----------------------------------|----------------|
| | Meat products | Dairy products | Meat products | Dairy products |
| Relative Accuracy | 94,7 % | 97,2 % | 94,0 % | 92,4% |
| Relative Sensitivity | 87,5 % | 96,8 % | 87,5 % | 85,7% |
| Relative Specificity | 100,0 % | 97,6 % | 100,0 % | 100 % |
| Selectivity | 100% | 100% | 100% | 100% |
| Detection Limits | 5- 10 UFC/ g | 50- 100 UFC/ g | 5- 10 UFC/ g | 50- 100 UFC/ g |

FACTORS THAT MAY INVALIDATE THE RESULTS

- Inaccurate standardisation of the inoculum.
- Material examined unsuitable.
- Use of expired systems.
- Temperature and time of incubation other than those recommended.

LIMITS AND WARNINGS

The presence of *Salmonella* spp., *Listeria* spp., *Escherichia coli*, *Staphylococcus aureus* and *Bacillus cereus* must be confirmed using appropriate culture media and biochemical and serological tests. For definitive identification of these micro-organisms, it is necessary to perform supplementary confirmation tests.

PRECAUTIONS

The product **FOOD SYSTEM** cannot be classified as hazardous under current legislation, nor does it contain harmful substances in concentrations $\geq 1\%$; it therefore does not require a Safety Data Sheet to be available.

FOOD SYSTEM must be used in the laboratory by properly trained personnel using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes of temperature. In such conditions the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.


DISPOSAL OF USED MATERIAL

After use, **FOOD SYSTEM** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

| Product | Code | Package |
|--------------------|-------|---------|
| FOOD SYSTEM | 71680 | 40 test |

TABLE OF SYMBOLS

| | | | | |
|---|---|--|---|--|
|  Batch code |  Do not reuse |  Manufacturer |  Contains sufficient for <n> tests |  Temperature limitation |
|  Catalogue number |  Fragile, handle with care |  Use by |  Caution, consult accompanying documents | |

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1. Difco & BBL Manual, Manual of Microbiological Culture Media 2003.
2. **ISO 6579:2002**. Horizontal method for the detection of *Salmonella* spp.
3. **ISO 6888-2:1999**. Enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species). - Technique that uses Rabbit Plasma Fibrinogen (RPF) Agar Medium.
4. **ISO 7932:2004**. Horizontal method for the enumeration of presumptive *Bacillus cereus*.
5. **ISO 11290:2002**. Horizontal method for the detection and enumeration of *Listeria monocytogenes*.
6. **ISO 16140:2003**. Microbiology of food and animal feeding stuffs - Protocol for the validation of alternative methods.

Microbiology Products

LIOFILCHEM S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net